



Diamant Schleifer für die Zahnmedizin  
 Diamond instruments for the dentistry  
 Instruments diamantés pour le odontologie  
 Instrumentos diamantados para odontologia

Anwendung / Application / Application / Empleo

Einsatzgebiet Application area Domaine d'application Usu en sector	Arbeitsgang Working process Phase de travail Fase de trabajo	Körnung (30x) Grit (30x) Grain (30x) Grano (30x)	Instrumentenbeispiel Examples of instr. Exempl. d'instruments Ej. de los instrumentos	ISO-Ø	empf. Drehzahl min <sup>-1</sup> Recommended speed Vit. de rotation rec. Veloc. recomendada
Schmelz, Dentin Enamel, Dentine Email, Dentine Esmalte, Dentina	Zügiges Vorschleifen faster pregrinding Dégrossissage rapide Rectificado inicial rápido		5831/014 	012-023	225.000 - 60.000
Schmelz, Dentin Enamel, Dentine Email, Dentine Esmalte, Dentina	Vorschleifen Pregrinding Dégrossissage Rectificado inicial		CD72G/014  6831/014 		
Schmelz, Dentin Enamel, Dentine Email, Dentine Esmalte, Dentina			CD34/012 	010-023	225.000 - 60.000
Schmelz, Dentin, Füllungen Enamel, Dentine, Fillings Email, Dentine, Obturations Esmalte, Dentina, Obturaciones	Formschleifen Shape grinding Façonnage Rectificado de forma		CD80/014 		
Schmelz, Dentin Enamel, Dentine Email, Dentine Esmalte, Dentina			831/012 	008-055	225.000 - 35.000
Schmelz, Dentin, Füllungen Enamel, Dentine, Fillings Email, Dentine, Obturations Esmalte, Dentina, Obturaciones			859/010 		
Schmelz, Dentin, Füllungen Enamel, Dentine, Fillings Email, Dentine, Obturations Esmalte, Dentina, Obturaciones	Feinschleifen Fine grinding Abrasion de précision Rectificado fino		CD72F/014  8801/023 	007-031	100.000 - 20.000
Wurzelszement Root cementum Ciment radiculaire Cemento de raíces	Glätten Polishing Lissage Alisado		8889/010 		
Schmelz, Dentin, Füllungen Enamel, Dentine, Fillings Email, Dentine, Obturations Esmalte, Dentina, Obturaciones			862EF/012 		
Schmelz, Dentin (für Ceramic Inlays) Enamel, Dentine (for Ceramic Inlays) Email, Dentine (pour Inlays en céramique) Esmalte, Dentina (para ceramic Inlays)	Finieren Finishing Finissage Acabado		CI847EF/018 	007-029	100.000 - 20.000
Wurzelszement Root cementum Ciment radiculaire Cemento de raíces	Glätten Polishing Lissage Alisado		957EF/009 		
Zahnstein, Konkremente Calculus, concrements Plaque dentaire, dépôts Cálculo, concreciones	Prophylaxe Prophylaxis Prophylaxie Profilaxis		8407LEF/012 	012	8.000 - 12.000
Zirkonoxid Zircon oxide Oxyde de zirconium Óxido de zirconio	Trepanieren/Trennen Trepanation/cutting Trépaner/couper Trepanación/Separado		6880Z 	010-014	160.000 - 140.000
Schmelz, Dentin Enamel, Dentine Email, Dentine Esmalte, Dentina	minimal invasiv minimum contact invasive minimale invasiva minimal		MI806  MI8806 	007-014 007-014	150.000 - 120.000 120.000 - 100.000